Melpop Trolley





EPP top opening isothermal container with extractable handle. For the transport of hot, chilled and frozen meals.

- Expanded polypropylene (EPP) isothermal container.
- Capacity: 60 litres.
- Equipped with two castors and an extractable ergonomic handle for easy handling.
- Designed for use with eutectic plates, to extend the transport time of hot, chilled and frozen food.
- It allows to properly operate in HACCP environments.
- The materials used for the manufacturing are suitable for food contact over the whole container surface.
- The unit's design allows accurate cleaning and sanitising (EC Reg. 852/2004 –HACCP-) including using a dishwasher.
- Fully recyclable at the end of the operating life.
- Guaranteed for use at temperatures ranging between -40° and +120°C.

Model	Code	Colour	External dimensions mm	Internal dimensions mm	Capacity I	Weight Kg
Melpop Trolley	MPTR0001	black	420x420xh600	360x360xh460	60	1,60
Chilled eutectic plate	PEML0001	pink	360x270xh39			2,9
Frozen eutectic plate	PEML0002	blue	360x270xh39			2,9
USB Data logger	AFG002		80x33xh14			



Via Savigliano 34 - 12030 Monasterolo di Savigliano (CN) Tel. +39 0172812600 - Fax +39 017285991 - info@melform.com GPS: Lat 44.681363 - Long 7.624809 <u>www.melform.com</u>